

K-COM KNEE BRACE

by ORTEMA

Carbon Fibre and Titanium

Partially flexible thigh and lower leg shells. Only 1,5 mm thick shells in sandwich construction with the ability to adapt to the changing muscular relief.

Optimal anatomical fit
Gives optimal wearing comfort.

Torsion-resistant, non-slipping carbon fibre construction. Stabilizes the joint where the knee needs it most.

Extremely light combined with stability. Maximal stabilization of the knee with a minimum of weight.

Polycentric joints made of titanium
Developed to meet biomechanical criteria-flat, light and extremely stable.

Flexion and extension limitation
Flexible adjustable.

Counter-rotating Velcro® strap system. Enabling the K-COM knee brace a non-slipping fit without adding to unnecessary thickness.

Test winner
The K-COM was awarded first place comparing 14 of the most commonly used knee braces.

The individual knee brace which one does not feel in everyday life and in sports.

E-Mail: servicefertigung@ortema.com

Knee // Shoulder // Elbow // Hand // Face

100 % individual

Manufacturing a K-COM knee brace according to your needs and specifications at our establishment for the areas knee, shoulder, elbow, hand and face.

Extreme thin carbon fibre material gives optimal wearing comfort with a minimum of weight and maximal stability.



ORTEMA GmbH SERVICEPRODUCTION:

Kurt-Lindemann-Weg 10 • 71706 Markgröningen • Phone +49 (0) 7145 - 91 53 800
Fax +49 (0) 7145 - 91 53 980 • E-Mail: info@ortema.de

ORTEMA Branch:

Pforzheim Tel. +49 (0) 7231 - 13 96 667
Ludwigsburg Tel. +49 (0) 7141 - 99 68 720
Neckarsulm Tel. +49 (0) 7132 - 38 32 866

ORTEMA Bases:

NORTH: ORTHO.FAKTUM Orthopädie-Technik, 21376 Salzhäusen
Phone +49 (0) 4172 - 40 89 944

EAST: Kröger's Centrum für Technische Orthopädie GmbH, 03238 Massen-Niederlausitz
Phone +49 (0) 3531 - 79 90 73 93

WEST: edirtarenaWest powered by KTM, 53619 Rheinbreitbach
Termine: Tel. +49 (0) 7145 - 91 53 890

www.ortema.de



Supply *your* Patients with the Proven ORTEMA Technology



Orthopädie-Technik & Sanitätshaus

Our Service for Established Doctors and Medical Establishments:

Supply *your* Patients with the Proven ORTEMA Technology

Now with the
3D Scan technology



Let us manufacture a K-COM knee brace according to your needs and specifications at our establishment. We'll happy work together with your orthopedic technician or your establishment. For this, we only need a plaster cast or the scan data of the patient's leg. Within ten workdays (plus shipping time), your orthopedic technician will receive the test-ready K-COM knee brace. This allows you to offer your patients the proven ORTEMA K-COM knee brace care and support at the highest level. Are you interested in this type of service?

Then feel free to call us, send us a fax or write us a mail:

Phone: + 49 (0) 7145 - 91 53 800
 Fax: + 49 (0) 7145 - 91 53 980
 E-Mail: servicefertigung@ortema.com

Benefit from the advantages of the **ORTEMA SERVICEPRODUCTION** for the areas knee, shoulder, elbow, hand and face.

- ▶ Save time and money
- ▶ Bridge staff shortages
- ▶ Open up new business fields
- ▶ Use our know-how for *your* patients



IN 5 STEPS TO THE SUPPLY WITH AN INDIVIDUAL ORTHOSIS OF ORTEMA